

PATIENT SPECIFIC (DISPOSABLE) APPLICATION INSTRUCTIONS

FULLBACK AND STANDARD SLINGS - FITTING FROM A SEATED POSITION



Apply and push the Sling down - as above. Ensure that the base of the Sling is about an Inch further down and below the base of the spine.



Pull the legs of the sling through, under the thighs, and at a right angle with the long bone of the thighs. Check that the legs are not twisted.

FULLBACK AND STANDARD SLINGS - FITTING FROM A SEATED POSITION (Continued)

Cross the leg displacement straps as shown. Check that the sling is free from twists, and is smooth under the patient. Move the hoist into position and attach the leg displacement straps first.



Note for a 2 point spreader bar, it is particularly important that the leg displacement straps are attached first. Raise the hoist and then attach the shoulder displacement straps. The hoist may move towards the patient. Ensure that the patient's arms and elbows are placed within the confines of the sling as shown in the drawing.

www.gomax.com

FULLBACK AND STANDARD SLINGS - ALTERNATIVE SEATED POSITION

For closed leg positions, the legs of the sling may be crossed under both the legs of the patient.



Note: The longer the shoulder displacement strap, the more reclined the patient will be. For a more upright seated position, shorten the shoulder displacement strap, and increase the leg displacement strap.